

PART 7 DETECTION OF MALIC ACID

PRODUCT CODE	EQUIPMENT REQUIRED	
LP TLCJARSTSG - 250	* SCREW TOP, SQUARE JAR – 250ML	
LP TLCTWEZS - 125	* TWEEZERS (STAINLESS STEEL) – 125MM	
LP TLCSPCLPL	* SPRING CLIP - LARGE	
LP MSCYLG - 10	* MEASURING CYLINDER, GRAD.GLASS – 10ML	
LP FFUNG - 40	* FILTER FUNNEL, GLASS – 40MM	
LP BEKSQG - 100	BEAKER, SQUAT, GRAD.GLASS – 100ML (4 x 100ML BEAKERS REQ'D @ \$5.45 ea)	
LP TLCSHT – 2020S	THIN LAYER CHROMATOGRAPHY (TLC) SHEET (MACHEREY – NAGEL) *200 x 200mm TLC SHT (1 x TLC SHEETS REQ'D)	
LP TLC CAPTB - 100	* TLC CAPILLARY TUBES – 100/ VIAL	
LP TLC BLLDPE - HB	* BLACK LEAD PENCIL, HB	
LP TLC RULCPER	* RULER, CLEAR PERSPEX – 30CM	
LP TLC WHTPPER – A4	* WHITE PAPER – A4 PAPER	
	REAGENTS	
LC TLCSTDM - 25	* TLC STANDARD MALIC ACID (25ML) – 2G/L	
LC TLCSTDL - 25	* TLC STANDARD LACTIC ACID (25ML) – 2G/L	
LC TLCSLMX - 100	* TLC SOLVENT MIXTURE – 100ML	
LC TLCINDIC - 100	* TLC INDICATOR – 100ML	